



Refined Strength | Advanced Aesthetics | Easy Debonding

# symetri™ CLEAR



†Data on file.

## Proprietary Base Design and Patented Debonding Plier for One-piece Removal



Ormco's patented laser etched pattern on the surface of the pad is precisely controlled to provide consistency and reliability in bonding and debonding.

Minimal force debonding with our patented debonding plier is designed to remove the bracket in one piece, without fracture. One squeeze. One click.

**Ormco™**

# symetri™ CLEAR

Part of Ormco's Straight-Wire™ system and designed with optimized in/out dimensions to express the purest McLaughlin, Bennett, Trevisi\* Rx.

Unique to Ormco - the optional slim in/out lateral aligns the incisal edges of the centrals and laterals.

\*Does not imply endorsement. Computer simulated data on file.

tooth	torque	ang	rot	bracket type	part numbers			
MAXILLARY					.018 L	.018 R	.022 L	.022 R
Central	+17°	4°	0°	twin	418-1111	418-1110	419-1111	419-1110
Lateral	+10°	8°	0°	twin	418-1211	418-1210	419-1211	419-1210
Lateral - Optional Slim In/Out	+10°	8°	0°	twin	418-1213	418-1212	419-1213	419-1212
Cuspid	0°	8°	0°	twin	418-1311	418-1310	419-1311	419-1310
Cuspid - Hook	0°	8°	0°	twin	418-2311	418-2310	419-2311	419-2310
Bicuspid	-7°	0°	0°	twin	418-1510		419-1510	
Bicuspid - Hook	-7°	0°	0°	twin	418-2511	418-2510	419-2511	419-2510
MANDIBULAR					.018 L	.018 R	.022 L	.022 R
Anterior	-6°	0°	0°	twin	418-1010		419-1010	
Cuspid	0°	3°	2°	twin	418-1411	418-1410	419-1411	419-1410
Cuspid - Hook	0°	3°	2°	twin	418-2411	418-2410	419-2411	419-2410
1st Bicuspid	-12°	0°	0°	twin	418-1610		419-1610	
1st Bicuspid - Hook	-12°	2°	0°	twin	418-2611	418-2610	419-2611	419-2610
2nd Bicuspid	-17°	0°	0°	twin	418-1710		419-1710	
2nd Bicuspid - Hook	-17°	2°	0°	twin	418-2711	418-2710	419-2711	419-2710

Single Patient Kits				.018 L	.018 R	.022 L	.022 R
Symetri Clear Upper/Lower 5-5 With Hooks 3-4-5				twin	746-2101		746-2201
Symetri Clear Upper/Lower 3-3 With Hooks on 3				twin	746-2102		746-2202
Symetri Clear Upper 3-3 With Hooks on 3				twin	746-2103		746-2203
Symetri Clear Upper 4-4 With Hooks 3-4				twin	746-2104		746-2204
Symetri Clear Upper 5-5 With Hooks on 3				twin	746-2105		746-2205
Symetri Clear Upper 5-5 With Hooks 3-4-5				twin	746-2106		746-2206



**Symetri Clear Debonding Tool**

866-4025

**Symetri Clear Typodont**

717-0021

Round and Smooth for  
Comfort and Aesthetics

